

EKR-Set

Removing of crowns, popular solutions of acurata




Selection of different toothings and head sizes for all challenges in your practice. The one-piece construction significantly reduces the risk of breakage.

030 314 EKR SET

Set includes 2 articles each:

		L	Ø mm
	500 314 139 015 012	4,2	1,2
	500 314 140 MK 014	5,0	1,4
	500 314 140 RR 012	5,0	1,2
	500 314 140 SCC 016	5,0	1,6
	800 314 473 528 016	8,0	1,6



Crown cutter gold

For gold and precious metals.

MK / RR crown cutter

Our safe problem solvers for crowns of non-precious alloys, metal crown with veneer, reduced alloys and precious metals.

SCC multifunctional separator

Quick and efficient separation of metal crowns and amalgam fillings.

LARRY BLACK

Our specialist for high performance ceramics like zirconium oxide or lithium disilicate.

Tip: tangential working and slight sawing movement facilitates the separation of ceramics