

COMP-Set







Polish of composite fillings, popular solutions of acurata



Exclusively for your practice: flames, cups, Turbo Spin and Turbo Shine - the ideal shape for all points of use

803 204 COMP SET

Set includes one pieces of each:

		L	Ø mm
	803 204 030 505 060	7,0	6,0
	803 204 243 505 040	10,0	4,0
	803 204 243 505 030	6,0	3,0
	803 204 030 495 060	7,0	6,0
	803 204 243 495 040	10,0	4,0
	803 204 243 495 030	6,0	3,0

		L	Ø mm
	803 204 011 505 055	8,0	5,5
	803 204 377 505 100	1,9	10,0
	803 204 011 495 055	8,0	5,5
	803 204 377 495 100	1,9	10,0

recommended speed 5.000 – 8.000 min⁻¹

The perfect polish in just 2 steps!

Pinky Standard/ Pinky Highshine

Works on all compomers and composites.
Doesn't smear and is very efficient.

Turbo Spin

By the brush-like characteristics of the polishing body also deep fissures are reached – multiple use and extremely durable.

Turbo Shine

Interspersed with diamond grain, adapts flexibly to nearly all surfaces and allows easy working without contact pressure.



Limited Edition