

Aligner polishers

Ideal for rounding off and smoothing the edge areas of thermoforming splints

- Soft elastomer fiber polisher for splint polishing
- With low heat generation and long service life thanks to impregnated special fiber
- Ideal for rounding off and smoothing the edge areas of thermoforming splints for aligner therapy, sports mouthguards, snoring splints and soft-retaining relinings in full dentures
- Rounding off and smoothing temporary restorations
- Smoothing 3D printed splints
- Optimum speed: 5.000 – 8.000 min⁻¹



014 900 372 SFT 250

Ø mm / L 25,0 / 4,0

packaging unit = 25 pieces



NEW

Try our Softies in the starter set

acurata mandrell 8 mm



330 104 604 391 080

Ø mm 8,0

packaging unit = 10 pieces

Starter set

014 900 372 SFT SET contains:

- 8x 014 900 372 SFT 250
- 1x 330 104 604 391 080

